



Recombinant Human Uncharacterized protein C9orf117 (C9orf117)

Product Code	CSB-EP686386HU-B
Abbreviation	C9orf117
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JU67
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAPKKSVSKA GKELEVKKKG GKKEPVAVE PPLAKEMKEF YHIQIRDLED RLARYQRKWD ELAVQEKMFR QEFEQLANNK KEIVAFKRT LNQQVDEITD LNEQLQNLQL AKEMEKDAFE AQLAQVRHEF QETKDQLTTE NIILGGKLA LEEFRLQKEE VTDKFTLLEE QVRKQENEFY DYAYNLEKKS VLDKDRLRKE IIQRVNLVAN EFHKVTTNRM WETTKRAIKE NNGITLQMAR VSQQGMKLLQ ENEQLKGRQN NLCKQLELLE NTQKVMARHK RGHQKIILML TKKCQEQQQD TKEAEELRLL LSQLEQRSQ LQVDNQALKS QRDQLSLQLE QQQVDLQRLQ QELANEQKVR ASLEAALVQA TSFLQNILQM HRDEEDSDVD VTFQPWHKEM LQQLLVMLSS TVATRPQKAA CPHQESQSHG PPKESRPSIQ LPRTGSLLPQ LSDITPYQPG DLGLVPRQVH IPPNPQDLRL LSYITRVGTF RAHSSPEMRA PGSLKRLEKF SLPEVPLRPK
Source	E.coli
Target Names	CFAP157
Protein Names	Recommended name: Uncharacterized protein C9orf117
Expression Region	1-520
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.