



Recombinant Human Migration and invasion-inhibitory protein (MIIP)

Product Code	CSB-BP686389HU
Abbreviation	MIIP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JXC2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVEAEELAQL RLLNLELLRQ LWVGQDAVRR SVARAASESS LESSSSYNSE TPSTPETSST SLSTSCPRGR SSVWGPPDAC RGDLRDVAR S GVASLPPAKC QHQESLGRPR PHSAPSLGTS SLRDPEPSGR LGDPGPQEAQ TPRSILAQQS KLSKPRVTFS EESAVPKRSW RLRPYLGYDW IAGSLDTSSS ITSQPEAFFS KLQEFRETNK EECICSHPEP QLPGLRESSG SGVEEDHECV YCYRVNRRLF PVPVDPGTPC RLCRTPRDQQ GPGTLAQPAH VRVSIPLSIL EPPHRYHIHR RKSFDASDTL ALPRHCLLGW DIFPPKSEKS SAPRNLDLWS SVSAEAQHQK LSGTSSPFHP ASPMQMLPPT PTWSVPQVPR PHVPRQKP
Source	Baculovirus
Target Names	MIIP
Protein Names	Recommended name: Migration and invasion-inhibitory protein Alternative name(s): IGFBP2-binding protein Invasion-inhibitory protein 45 Short name= Ilp45
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.