



Recombinant Arabidopsis thaliana 60S ribosomal protein L8-2 (RPL8B)

Product Code	CSB-BP686901DOA
Abbreviation	RPL8B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4PSL7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGRVIRAQRK GAAGSVFKSH THHRKGPAKF RSLDYGERNG YLKGLVTEII HDPGRGAPLA RVAFRHPFRY MKQKELFVAA EGMYTGGYLY CGKKANLMVG NVLPLGSIPE GAVICNVELH VGDRGALARA SGDYAIVIAH NPESNTTRVK LPSGSKKILP SACRAMIGQV AGGGRTEKPL LKAGNAYHKY KAKRNCWPVV RGVAMNPVEH PHGGGNHQHI GHASTVRRDK SAGAKVGQIA ARRTGRRRGA AATLAAKADY
Source	Baculovirus
Target Names	RPL8B
Protein Names	Recommended name: 60S ribosomal protein L8-2
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.