



Recombinant *Macaca fascicularis* STE20-related kinase adapter protein alpha (STRADA)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-BP687020MOV |
| Abbreviation | STRADA |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q4R6X5 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSFLTNDASS ESIAFSKQE VMSSFLPEGG CYELLTVIGK GFEDLMTVNL ARYKPTGEYV TVRRINLGAC SNEMVTFLQG ELHVSKLFNH PNIVPYRATF IADNELWVVT SFMAYGSAKD LICHTFMDGM NELAIAYILQ GVLKALDYIH HMGYVHRSVK ASHILISLDG KVYLSGLRSN LSMISHGQRQ RVVHDFPKYS VKVLPWLSPE VLQQNLQGYD AKSDIYSVGI TACELANGHV PFKDMPATQM LLEKLN GTVP CLLDTSTIPA EELTMSPSRS VANSGLSDSL TTSTPRPSNG DSPSHPYHRT FSPYFHHFVE QCLQPNPDAR PSASTLLNHS FFKQIKRRAS EALPELLRPV TPITNFESSQ SQDHSGIFGL VTNLEEELEVD DWEF |
| Source | Baculovirus |
| Target Names | STRADA |
| Protein Names | Recommended name: STE20-related kinase adapter protein alpha Short name=STRAD alpha Alternative name(s): STE20-related adapter protein |
| Expression Region | 1-394 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |