



Recombinant Human Synaptotagmin-like protein 3 (SYTL3)

Product Code	CSB-EP687190HU-B
Abbreviation	SYTL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4VX76
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAQEIDLSAL KELERAILQ VLYRDQAVQN TEEERTRKLIK THLQHLRWKG AKNTDWEHKE KCCARCQQVL GFLLRGAVC RGCShRVCAQ CRVFLRGTHA WKCTVCFEDR NVKIKTGEWF YEERAKKFPT GGKHETVGGQ LLQSYQKLSK ISVVPPTPPP VSESQCSRSP GRLQEFQFR GFNKSVENLF LSLATHVKKL SKSQNDMTSE KHELLATGPRQ CVGQTERRSQ SDTAVNVVTR KVSAPDILKP LNQEDPKCST NPILKQQNLP SSPAPSTIFS GGFRHGSLIS IDSTCTEMGN FDNANVTGEI EFAIHYCFKT HSLEICIKAC KNLAYGEEKK KKCNPYVKTY LLPDRSSQ GK RKTGVQRNTV DPTFQETLKY QVAPAQLVTR QLQVSVVHLG TLARRVFLGE VIIPLATWDF EDSTTQSFRW HPLRAKAEKY EDSVPQSNGE LTVRAKLVLP SRPRKLQEAQ EGTDQPSLHG QLCLVVLGAK NLPVRPDGTL NSFVKGCLTL PDQQKLRLKS PVLKQACPQ WKHSFVFSGV TPAQLRQSSL ELTVWDQALF GMNDRLLGGT RLGSKGDTAV GGDACSLSKL QWQKVLSSPN LWTDMTLVLH
Source	E.coli
Target Names	SYTL3
Protein Names	Recommended name: Synaptotagmin-like protein 3 Alternative name(s): Exophilin-6
Expression Region	1-610
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.