



# Recombinant Danio rerio Adenylate kinase 8 (ak8)

<b>Product Code</b>	CSB-MP687353DIL
<b>Abbreviation</b>	ak8
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q502L7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDSAERPLRI PPEMVIYAEK HDIFHLIQTL VRNLMVDKPD DPIQYLINFL ERDNVDVPRI ILLGPPASGK KTVAKKLCEH TQAIHITFCD ILKDDSDLTR AAQSYDKKQ NVPKDLWIQL IQQRLSKPDC VQRGWVLEAI PKTQDEALCL QEAGITPDHV VMLEAPDVVL IERSSGKRID PVTGDVYHVT FMWPESEEEVA QRLETPRTL M PVQVMSQKLQ KYHTEAHALK RTYHNCLKII NADQPHVDVF SQVLNFICTP RHSPSPYTPR ILLFGPPGAG RNLQAKLIAQ KYGIINICCG ELLKAUSADE SHMGELIKPY LESEQQVPSS IVLRLVTERL SRMDCTARGW VLHGFP RDVE EAELLHKS NF TPNRVFFLEI TDDIAIERLA LRAVDPVTGE WYHSVYKPPP GPEVQARLRF NPQHSEAQLL MRLKEYWSQT VRLQAFYPQA VYINADQDPH TVFESLESRL VGQLPKVLSN QQTNEN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ak8
<b>Protein Names</b>	Recommended name: Adenylate kinase 8 EC= 2.7.4.3
<b>Expression Region</b>	1-486
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.