



Recombinant Protein IncC (incC)

Product Code	CSB-BP687428ENL
Abbreviation	incC
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q52312
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	≥85% (SDS-PAGE)
Sequence	<p>MGAIHEETAN RSPIDGHHQ AGDRAADHRH SARRAGRWPA PGGVCHFAGA DQGGGVAGGQ SRMGSSRGTA PRGLRASDGG TAGASSVHRQ EVGSRRQEET GTQVMKTLVT AIQKGGQGKT FATCHLAFDF LERGLRVAVI DLDTQGNASF TLSAYQSGYL ASQLFTGDTD DLRYWFGKRE GESLALIAAD ANLANLDKME LSQAAAALRA SVAALGEFFD VCLIDTAPSL GVAMTAAVLT ADYMLSPIEM EAYSLQGMKK MVAVISNLRK QNPKLRFLGM VPNKVDARKP RHVNNLATLQ QAYPQLILPF SIGARDSIAE ALGEQMPVWK IKKTAARKAT QEVRALADYV FTKMEIAQ</p>
Source	Baculovirus
Target Names	incC
Protein Names	Recommended name: Protein IncC
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>