



Recombinant Human Raftlin-2 (RFTN2)

Product Code	CSB-EP687456HU
Abbreviation	RFTN2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q52LD8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GCGLRKLED PDDSSPGKIF STLKRPQVET KTEFAYEYVL LDFTLQASSN PEVIKINSIL DIVTKVENYY LKGYIVGAIH PVIQPVGQRK HLPASYLYRV VLLRLKLSPK NSAAPSGQRR PRLVIEECPL TSEAQTNDAA KELIEKINVA AKRGMKFVGF ISQHYSKSKF CNGTNHDGDI ESMLHVRHGS DENCRSWNEG TLGQSSSEG IEEELHHESG QYQMEQNGSP TSSKSRKGEA SDNKLYTVFN AFDDDSTSWA YQEGILSMKV TRKGSVISTL DADWLELTF YYKQGLSLID SFVFWETSKG EHLPKSLEGF FIYEEEGSGV PGSSRKGNDV IVVEQWTVIE GCEIKTDYGP LLHTLAEFGW LLTSVLPTPV LRHDSEGNLA TKQIVFLQRP VMWNSAAQTP DKKASRHIKG EDKNKATSR IGLDTTSSQP AESRHLPEEC RLSPSRECWT KEGRLAQHNS FSGFSSSDNV LREDDGQFD QEDGVTQVTC M
Source	E.coli
Target Names	RFTN2
Protein Names	Recommended name: Raftlin-2 Alternative name(s): Raft-linking protein 2
Expression Region	2-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.