



Recombinant Human Raftlin-2 (RFTN2)

Product Code	CSB-MP687456HU
Abbreviation	RFTN2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q52LD8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	GCGLRKLED PDDSSPGKIF STLRPQVET KTEFAYEYVL LDFTLQASSN PEVIKINSIL DIVTKVENYY LKGIVGAIH PVIQPVGQRK HLPASYLYRV VLLRLKLSPK NSAAPSGQRR PRLVIEECPL TSEAQTNDAA KELIEKINVA AKRGMKFVGF ISQHYSKSKF CNGTNHDGDI ESMLHVRHGS DENCRSWNEG TLGQSSSESG IEEELHHESG QYQMEQNGSP TSSKSRKGEA SDNKLYTVFN AFDDDSTSWA YQEGILSMKV TRKGSVISTL DADWLELTF YYKQGLSLID SFVFWETSKG EHLPKSLEGF FIYEEEGSGV PGSSRKGNDV IVVEQWTVIE GCEIKTDYGP LLHTLAEFGW LLTSVLPVTPV LRHDSEGNLA TKQIVFLQRP VMWNSAAQTP DKKASRHIKG EDKNKATSRN IGLDTTSSQP AESRHLPEEC RLSPSRECWT KEGRLAQHNS FSGFSSSDNV LRELDGQFD QEDGVTQVTC M
Source	Mammalian cell
Target Names	RFTN2
Protein Names	Recommended name: Raftlin-2 Alternative name(s): Raft-linking protein 2
Expression Region	2-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.