



Recombinant Dictyostelium discoideum DnaJ homolog subfamily A member 1 homolog (dnaja1)

Product Code	CSB-YP687550DKK
Abbreviation	dnaja1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q54ED3
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MVKEKEYYER LGVKPDCTED ELKKAYRKMA VKYHPDKNQG PGKDAAEAKF KDISEAYEVL SDPEKRKMYD SYGSEGMKES GFHASSAEDL FSHFFGAGGG GGGFSFGGGG GDDFGGFSFG NMGGMGGMGG MGGGHKKRRK GEDIAHEMNR SLEELYNGKL VKISISRDEV CKTCKGSGSN KPGVTTTCPT CNGSRYVFQK KQVGPGMIQQ VQTACHTCHG TGEKIKEEDK CKECKGKRVI QGKKIVQFQV EKGTRDGERI MLQGQGSEYP GVPPGDVIIT IREKPNVNFK RNGDNLIYTK RLKLLDSIAG SQFIINTLDQ RKLWVNHEKG DIIKQGMRY IENEGMPIKG TSRKGLIIA FDIEYPSNLT NDDIEKLSKI LPKAATPSVS KSDCKSVGLS KVNFTNEQS SHGGAGGAYQ QHGGAYGHQK QQQQGFNPAD FGAQFGGGGP QQAQQC
Source	Yeast
Target Names	dnaja1
Protein Names	Recommended name: DnaJ homolog subfamily A member 1 homolog
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.