



Recombinant Dictyostelium discoideum Protein arginine N-methyltransferase 1 (prmt1)

Product Code	CSB-EP687551DKK-B
Abbreviation	prmt1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q54EF2
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MTETNKNVDA LEKINQMSSA DYYFDSYSHF GIHEEMLKDE VRTLAYRRAI INNRKLFEGK VVLDVGCCTG ILCMFAAQAG AKMVIGVDNS EMLPIAQKII TANNFDKTIT LIKGKMEEVV LPVDKVDIII SEWMGYFMLY EGMLDVTLYA RDKYLVPGGV ILPDKASLYI TAIEDQDYKE EKINYWNNVY GFDMSCIREI ALKEPLVDVV QPNMIVTND CILTVDIMTI TKDELKFRSD FKLKALRDDL IHAFVVYFDI EFSKGDKPVC FSTGPKAKYT HWKQSIMYFE DHIKIQQGEI ITGTMDCAPF DKNQRDLKIK LDFNFAGELM KSSSSLEYHM R
Source	E.coli
Target Names	prmt1
Protein Names	Recommended name: Protein arginine N-methyltransferase 1 EC= 2.1.1.- Alternative name(s): Histone-arginine N-methyltransferase PRMT1 EC= 2.1.1.125
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.