



Recombinant Dictyostelium discoideum Checkpoint protein hus1 homolog (hus1)

Product Code	CSB-MP687659DKK
Abbreviation	hus1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q54NC0
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MKFKAKISKA SILLKTVQNI MKIYNECICH ITPDKLRFII QSDFNDGMQV FCEIQRSLIF ENFTIESLSD NEIQFQLDLE SLRRVLQSAT SNGVSDIFTN LTKVRGGPVL HFIKSGTSG TVIFQDIPIV LLTAQQMAQI NEPTLPDPLV HILLPNLKNL QKVIDKMKNI SDCLKIMIAM NNRLSFEVET SSGSISTFYS GLDHPQFGDH VLSPDLQATV CVDIKKFAKV LHIHQLKPSE VVLCLYERSI IVHVVLTDIM ITYYLPVLIK N
Source	Mammalian cell
Target Names	hus1
Protein Names	Recommended name: Checkpoint protein hus1 homolog
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.