



Recombinant Rat Protein BTG2 (Btg2)

| | |
|--------------------------|---|
| Product Code | CSB-MP002864RA |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P27049 |
| Product Type | Recombinant Protein |
| Immunogen Species | Rattus norvegicus (Rat) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSHGKRTDML PEIAAAVGFL TLLRTRGCV SEQRLKVFSR ALQDALTDHY KHHWFPEKPS KSGSYRCIRI NHKMDPIISK VASQIGLSQP QLHQLLPSEL TLWVDPYEVS YRIGEDGSIC VLYEEAPVAT SYGLLTCKNQ MMLGRSSPSK NYVMTVSS |
| Source | Mammalian cell |
| Target Names | Btg2 |
| Protein Names | Recommended name: Protein BTG2 Alternative name(s): BTG family member 2 NGF-inducible anti-proliferative protein PC3 |
| Expression Region | 1-158 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Target Details | This protein is a member of the BTG/Tob family. This family has structurally related proteins that appear to have antiproliferative properties. This encoded protein is involved in the regulation of the G1/S transition of the cell cycle. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |