



Recombinant Dictyostelium discoideum Serine/threonine-protein kinase 19 homolog (STK19)

Product Code	CSB-YP687756DKK
Abbreviation	STK19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q54WC0
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MNIKRNQNNNS IDNNLSIDSN PSPIKKVKNH LDFSKSYNED DEPSNDLDLS FNTSNISNLS NINGEEDNDD DDRENEPKDI DNPNPSTTT TATSTNSTLV KKTNNFLDFN FAFEPNKEED DEEEDLNISL TTEDILLSQI ELLKAEWDYI TYNKSPPLVL QSVYHQSNKD SERSIQILKK KGTIRLFQIA TSLNDYCIII TTDYIKNIEL LKIENTTIPN KIKPETLTNN MKNSSPIKSA PIDTILLKRS EILDLFIKKL IPNFNEVSIT RSKLFQLLSI VNDDHHQQEN IITHLVQCGL LLQKDDCTFT FSVPGAGGFF LNLMKGRKEI LSNIQRLKYK EILKKDLLKK KLKYSNVQMQ LLIKDLLGLN KIKIIPPTTQG ELIRCVQFDE L
Source	Yeast
Target Names	STK19
Protein Names	Recommended name: Serine/threonine-protein kinase 19 homolog EC= 2.7.11.1 Alternative name(s): STK19 homolog
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.