



Recombinant Dictyostelium discoideum Succinate dehydrogenase [ubiquinone] iron-sulfur subunit, mitochondrial (sdhB)

Product Code	CSB-YP687870DKK
Abbreviation	sdhB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q55CC2
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	TTTAKTE ASSSPQMKIK TAEKKDHFVN FQVYRYNEET TAKPYIQTYN INLKDCGPMV LDALLLIKNN IDPTLSFRRS CREGICGSCA MNLNGSNTLA CTKKITDCLS GDTVKVYPLP HMHVVRDLIP DLSHFYTQHK SIKPWLEPAV DVPRYNGKEI LQSKENRHKL DGLYECILCA CCSTSCPSYW WSEGGDGGYL GPAVLLQAYR WIADSRDSIQ KDRLAILED QMKVYKCHTI MNCTAVCPKG LNP GKSIANI KYLLAHN
Source	Yeast
Target Names	sdhB
Protein Names	Recommended name: Succinate dehydrogenase [ubiquinone] iron-sulfur subunit, mitochondrial EC= 1.3.5.1 Alternative name(s): Iron-sulfur subunit of complex II Short name= lp
Expression Region	24-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.