



Recombinant Dictyostelium discoideum Protein SUGT1 homolog (sugt1)

Product Code	CSB-MP687888DKK
Abbreviation	sugt1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q55ED0
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MEQLKEGNSY FVDEQYDEAL KCYDKACLEL SNNAEAFFKR SQCHSKLSNL KEALSDINTS IKLDSNNSKY YLKKGQLCFE LEEFDTALKT FEKGQSIDSE NSSFKTWIRK SKAEIQSNPT TTTTTTPTPT PTPTPAPQPV TTTTNTPIPI TTSNTTTTTN NNNNNNNNNN NNNNNNNTTT DSTTTKLPPI SSGNKVRHEW YQTETHVVL IFAKFVTASN SKINLTSKSV NISFPLANGS EFLFEMDLFD PIVDKDSTIH YYSTKVEIKM KKSRAIKWDT LEFTDKSGPV GLMDQISSSP AVSPYASKK DWDKLPNEPE EKLEGDQALN KIFRDIFSKG SEDQQRAMMK SFTESGGTVL STNWDEVGSK KVVGEPPKGL EFKQYEK
Source	Mammalian cell
Target Names	sugt1
Protein Names	Recommended name: Protein SUGT1 homolog
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.