



# Recombinant Mouse DEP domain-containing mTOR-interacting protein (Deptor)

<b>Product Code</b>	CSB-EP687983MO-B
<b>Abbreviation</b>	Deptor
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q570Y9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEEGSSGGSG SSSDSNAGGSG GVQQRELERM AEVLVTGEQL RLRLHEEKVI KDRRHHLKTY PNCFVAKELI DWLIEHKEAS DRETAIKLMQ KLADRGIIHH VCDEHKEFKD VKLFYRFRKD DGTFALDSEV KAFMRGQRLY EKLMSPETTL LQPREEEGVK YERTFMASEF LDWLVQEGEA TTRKEAEQLC HRLMDHGIIQ HVS NKHPFVD SNLLYQFRMN FRRRRRLMEL LNETSPSSQE THDSPFLCRK QSHDSRKSTS FMSVSPSKEI KIVSAVRRSS MSSCGSSGYF SSSPTLSSSP PVL CNPKSVL KRPVTSEELL TPGAPYARKT FTIVGDAVGW GFVVRGSKPC HIQAVDPSGP AAAAGMKVCQ FVSVNGLNV LNVDYRTVSN LILTGPRIV MEVMEELDC
<b>Source</b>	E.coli
<b>Target Names</b>	Deptor
<b>Protein Names</b>	Recommended name: DEP domain-containing mTOR-interacting protein Alternative name(s): DEP domain-containing protein 6
<b>Expression Region</b>	1-409
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.