



# Recombinant Rabbit Peptidyl-prolyl cis-trans isomerase FKBP4 (FKBP4)

<b>Product Code</b>	CSB-EP008700RB-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P27124
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Oryctolagus cuniculus (Rabbit)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTAEEMKAAE SGAQSAPLPL EGVDISPKQD EGVLKVIKRE GTGTETPMIG DRVHVHTGW LLDGTFDSS LDRKDKFSFD LGKGEVIKAW DIAVATMKVG ELCRITCKPE YAYGSAGSPP KIPP NATLVF EVELFEFKGE DLTDEEDGGI IRRIRTRGEG YARPNDGAIV EVALEGGYKD RLFQDRELRF EVGEGESLDL PCGLEKAIQR MEKGEHSILY LKPSYAFGNA GKEKFQIPPY AELKYEVLHK SFEKAKESWE MSSEKLEQS AIVKERGTVY FKEGKYKQAL LQYKKIVSWL EYESSFSSEE VQKAQALRLA SHLNLAMCHL KLQAFSAAVE SCNKALELDS NNEKGLFRRG EAHLAVNDFD LARADFQKVL QLYPSNKA AK AQLAVCQQRI RKQIAREKKL YANMFERLAE EENKAKAEVA AGDHPMDTEM KDERNDVAGS QSQVETEA
<b>Source</b>	E.coli
<b>Target Names</b>	FKBP4
<b>Protein Names</b>	Recommended name: Peptidyl-prolyl cis-trans isomerase FKBP4 Short name= PPIase FKBP4 EC= 5.2.1.8 Alternative name(s): 52 kDa FK506-binding protein Short name= 52 kDa FKBP Short name= FKBP-52 59 kDa immunophilin
<b>Expression Region</b>	1-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This protein is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. This encoded protein is a cis-trans prolyl isomerase that binds to the immunosuppressants FK506 and rapamycin. It has high structural and functional similarity to FK506-binding protein 1A (FKBP1A), but unlike FKBP1A, this protein does not have immunosuppressant activity when complexed with FK506. It interacts with interferon regulatory factor-4 and plays an important role in immunoregulatory gene expression in B and T lymphocytes. This encoded protein is known to associate with phytyl-CoA alpha-hydroxylase. It can also associate with two heat shock proteins (hsp90 and hsp70) and thus may play a role in the intracellular trafficking of hetero-oligomeric forms of the steroid hormone receptors. This protein correlates strongly with adeno-associated virus



type 2 vectors (AAV) resulting in a significant increase in AAV-mediated transgene expression in human cell lines. Thus this encoded protein is thought to have important implications for the optimal use of AAV vectors in human gene therapy. The human genome contains several non-transcribed pseudogenes similar to this gene.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.