



Recombinant Methanocaldococcus jannaschii UPF0252 protein MJ1282 (MJ1282)

Product Code	CSB-YP688299MRU
Abbreviation	MJ1282
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q58678
Product Type	Recombinant Protein
Immunogen Species	Methanocaldococcus jannaschii (strain ATCC 43067 / DSM 2661 / JAL-1 / JCM 10045 / NBRC 100440) (Methanococcus jannaschii)
Purity	>85% (SDS-PAGE)
Sequence	MRKIILIFFI FLILQNIYAY EKIYDVNHVC YYTLTSVELQ KFKVKTAEIIT PLNADNKTTI TKDTILVGNP EENPLTKKYI GFFKIKINKT FPGKNKGIIE KQIINGHTVI LLAGSDIQGT YASITFANL NDIPENPIIC ETSNKVNIYS TSLNSEYFRD FIEKNILTPK EIEKVKLSY KLKGDKKAT IENIAKWVAN NIKYDYDKCK KIESGKFSWD EYNTPSETVK TKGVCCLDYA TLTSALLND NIIPYMLDVA LYNTSSLKIS SYHASVAVKI DNTYFVIDQQ PYLIPINEYT AQT FEDNLR IASIVMFRVVK ERDGIKLIKE KEVPGIAIYG DLINLLEMRF NN
Source	Yeast
Target Names	MJ1282
Protein Names	Recommended name: UPF0252 protein MJ1282
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.