



Recombinant Bovine Sestrin-2 (SESN2)

Product Code	CSB-EP688343BO-B
Abbreviation	SESN2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q58CN8
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	MIVADSECRA ELKGYLPAG EEQRESRVRR GPRGPSAFIP VEEVLQEGAE SLEQHLGLEA LMSSGRVDNL AVVMGLHPDY FTSFWRLHCL LLHTDGPLAN SWRHYIAIMA AARHQCSYLV GSHMAEFLQT GGDPEWLLGL HRAPEKLRKL SEINKLLAHR PWLITKEHIQ ALLKTGEHSW SLAELIQALV LLTHCHSLAS FVFGCGILPE GDPEGSPAPQ APSPPEQST PPSRDSL NHS GGFEAARDVE ALMERMRQLQ ESSLQDEGAS QEEMESRFEL EKSESLLVTP SVDILEPFAN PDMLCFVEDP TFGYEDFTRR GTQAPPTFRA QDYTWEDHGY SLIQRLYPEG GQLLDEKFQA AYSLTYNTIA MHSGVDTSVL RRAIWNYIHC VFGIRYDDYD YGEVNQLLER NLKVYIKTVA CYPEKTTRRM YNHFWRHFRH SEKVHVNLLL LEARMQAALL YALRAITRYM T
Source	E.coli
Target Names	SESN2
Protein Names	Recommended name: Sestrin-2
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the sestrin family of PA26-related proteins. The encoded protein may function in the regulation of cell growth and survival. This protein may be involved in cellular response to different stress conditions.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.