



Recombinant Danio rerio Isthmin-2 (ism2)

Product Code	CSB-YP688371DIL
Abbreviation	ism2
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q58T08
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	FPVRH HKRSRNGVYV ENQVQNHLAD SHPHQRRWLQ HHSTGVLPLP EPEEESKPFV LDFKNLPDLA NADIGSQNPN IQVTIEVLDD PPMDVEMDLV KEWSNDWSTS SPSSTVEWLG GKCLFWPLFW GYTADDSGED GTGQAEDEED DYDYDSGEPI PSGLGKTDGD WTHNRYEEKE EEWSTWSPCS VTCGHGNQTR SRSCGDFCTS TESQSCDLVP CPDDWNSVGH VPFEMENGT EPYGTVDVGS EKWLNCKSEF LQRYLQQVFT ELPNCPCSY SHVSNNIVSL LDVGHERSFQ WRDASGPKER LDIYKPSARS CLRSGLSKNG TTLAAQHCCY DDNKQLITRG KGAGTPNLIS TEFSPELHFK VDVLPWILCK GDWSRFHAVR PPNNGLHCME NPQQDIFMNE LEEAREY
Source	Yeast
Target Names	ism2
Protein Names	Recommended name: Isthmin-2
Expression Region	26-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>