



Recombinant Candida albicans Pre-mRNA-splicing factor ISY1 (ISY1)

Product Code	CSB-BP688435CZD
Abbreviation	ISY1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q59R35
Product Type	Recombinant Protein
Immunogen Species	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSRNKEKAQS SLNRFYQSQS SAPINYHYHQ RPKNIHSITQ LSQAEGWRRS IIGEISKQLT EINDFEISDI KVRELNDQLN ELFNEKRRWE YHIRDKLHGN DYIRHNGNTK NDLINTGIRV NNSDKGKYR YFGRAKDLPE VKLMIEENKK NRKNHKRDDI DDSQSKRRR KVDQYKSTDN SIDMNYGIFY DEEEDYLDLDD SEAQLGDPLV VYEEKRTKEL NDAEIQQQDT ENLEQIRMSF NDIPTNEVVT KWLVNKRREQ LLSKFR
Source	Baculovirus
Target Names	ISY1
Protein Names	Recommended name: Pre-mRNA-splicing factor ISY1
Expression Region	1-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.