



# Recombinant Danio rerio Ectonucleotide pyrophosphatase/phosphodiesterase family member 6 (enpp6)

<b>Product Code</b>	CSB-YP688554DIL
<b>Abbreviation</b>	enpp6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5BKW7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	NRKLLVFL IDGFRHDYMD DLHNLPGFRE IVENGVKVDY LTPDFPSSLY PNYYSLMTGR HCEVHQMTGN YMWDTDTQKE FLIGTNPDSR LPMWWDGSEP LWVTMQKLGK KVYMYWPGC EVTILGVRPT FCEEYVYNPS EKNLTDSMEN ALNALKSSKA DMAGIYYEKI DVEGHHFGPR SPEIQRAIRS LDQAFQILNQ KIREKNMRDT INVVLFS DHG MTQLKWMEKI IELDNYINMS HIIKMMDRGP VVSLWPKQDK FEEIQNLST ADNMMVYKKH EIPDRFHYKN GQFVSTLTLV AEPGWFITEN KAKLPFWNNG TEAAGGWQHG WHGYDNEFVD MRGSFLAQGP DFKSNYRAGP IRTVDVYNVL CKTLGMNPLP NNGSWSRVEC MMRSS
<b>Source</b>	Yeast
<b>Target Names</b>	enpp6
<b>Protein Names</b>	Recommended name: Ectonucleotide pyrophosphatase/phosphodiesterase family member 6 Short name= E-NPP 6 Short name= NPP-6 EC= 3.1.-.-
<b>Expression Region</b>	23-415
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Extracellular domain
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.