



Recombinant Bovine Heme oxygenase 1 (HMOX1)

Product Code	CSB-MP688656BO
Abbreviation	HMOX1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5E9F2
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	MERPQPDSSM PQDLSEALKE ATKEVHTQAE NAEFMKNFQK GELTQEGFKL VMASLYHIYV ALEEEIERNK ENPVYTPLYF PEELHRRASL EQDMAFWYGP RWQEAIPYTQ ATKRYVQRLQ EVGRTEPELL VAHAYTRYLG DLSGGQVLKK IAQKALNLP SGEGLAFFTF PNIASATKFK QLYRSRMNTL EMTPEVRQRV LDEAKTAFLL NIQLFEELQG LLTQKAKDHD PLQAPELHRR AGSKVQDLAP TKASRGKPQP SVLSQAPLLR WVLTL SFLVA TVAVGLYAM
Source	Mammalian cell
Target Names	HMOX1
Protein Names	Recommended name: Heme oxygenase 1 Short name= HO-1 EC= 1.14.99.3
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	Heme oxygenase, an essential enzyme in heme catabolism, cleaves heme to form biliverdin, which is subsequently converted to bilirubin by biliverdin reductase, and carbon monoxide, a putative neurotransmitter. Heme oxygenase activity is induced by its substrate heme and by various nonheme substances. Heme oxygenase occurs as 2 isozymes, an inducible heme oxygenase-1 and a constitutive heme oxygenase-2. HMOX1 and HMOX2 belong to the heme oxygenase family.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.