



Recombinant Bovine Hsp90 co-chaperone Cdc37 (CDC37)

Product Code	CSB-EP688677BO-B
Abbreviation	CDC37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5EAC6
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MVDYSVWDHI EVSDDEDETH PNIDTASLFR WRHQARVERM EQFQKEKEEL DRGCRECKRK VAECQRKLKE LEVAEGEGGK AELERLQAEA QQLRKEERSW EQKLEEMRKK EKSM PWNVDT LSKDGF SKSM VNTKPEQAE ESEEVREQKH KTFVEKYEKQ IKHFGMLRRW DDSQKYLSDN VHLVCEETAN YLVIWCIDLE VEEKCALMEQ VAHQ TIVMQF ILELAKSLKV DPRACFRQFF TKIKTADRQY MEGFNDELEA FKDRVRGRAK LRIEKAMKEY EEEERKKRLG PGGLDPVEVY ESLPEELQKC FDAKDVQMLQ DAISKMDPTD AKYHMQRCID SGLWVPNSKS SEAKDGEEAG TGDPLLEATS KSGDEKDVSV
Source	E.coli
Target Names	CDC37
Protein Names	Recommended name: Hsp90 co-chaperone Cdc37 Alternative name(s): Hsp90 chaperone protein kinase-targeting subunit p50Cdc37
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is highly similar to Cdc 37, a cell division cycle control protein of <i>Sacchomyces cerevisiae</i> . This protein is a molecular chaperone with specific function in cell signal transduction. It has been shown to form complex with Hsp90 and a variety of protein kinases including CDK4, CDK6, SRC, RAF-1, MOK, as well as eIF2 alpha kinases. It is thought to play a critical role in directing Hsp90 to its target kinases.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.