



Recombinant Danio rerio RING1 and YY1-binding protein B (rybbp)

Product Code	CSB-MP688689DIL
Abbreviation	rybbp
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5EG55
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MGDKKSPTRP KRQAKPTADN GFWDCSVCTF RNSAEAFKCS ICDVRKGTST RKPRINSQLV AQQVAQQYAT PPPPKKEKKE KPERPEKDRA EEERPDINPP DEHPVEQRDK DKSEKEQPEK EKKDREKEII PAITKKPNSK KNRPKSDIHQ SPPSERNSIQ SGKSTTKTKN SHNSRPKLKN IDRSTAQQLA ITVGNVTVII TDFKEKTRTS STSSSTVTSS ASSEQQHQSS GSESMDKGSS RASTPKGDL LGHDESF
Source	Mammalian cell
Target Names	rybbp
Protein Names	Recommended name: RING1 and YY1-binding protein B Alternative name(s): Death effector domain-associated factor B Short name= DED-associated factor B
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.