



Recombinant *Xenopus laevis* Maspardin (spg21)

Product Code	CSB-EP688860XBE
Abbreviation	spg21
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5FVD6
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MGAIKISPDY NWFRSTVPLK RIIVDDDDSK VWSLYDAGPR SVRCPIIYLP PVSQTADVFF HQILALTGWG YRVIALQYPV YWDHLEFCDG FRKLLDDLHL DKVHLFGASL GGFLAQKFAE YTHKSPRVQS LILCNSFSDT SIFNQTWTAN SFWLMPAFML KKILLGNFSS GPVDPVMADA IDFMVDRLES LNQSELASRL TLNCQNAYVE PHKIRDYAVT IMDVFDQSAL SMDAKEEMYK LYPNARRAHL KTGGNFPYLC RSAEVS LYVQ IHLRQFHGTR YAAIDPAMVS EEELEVQKSF PQTTVEAEES
Source	E.coli
Target Names	spg21
Protein Names	Recommended name: Maspardin Alternative name(s): Spastic paraplegia 21 autosomal recessive Mast syndrome protein homolog
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.