



# Recombinant Rat Proline-rich protein 5 (Prr5)

<b>Product Code</b>	CSB-BP688861RA
<b>Abbreviation</b>	Prr5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5FVG6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRTLRLKFM SSPSLSDLGK REPGAAGADE RGTQQRRACA NATWNSIHNG VIAVFQRKGL PDQELFILNE GVRQLLKTEL GSFFTEYLQN QLLTKGMVIL RDKIRFYEGQ KLLDSLAEW DFFFSDVLP LQAIFYPVQG KEPSVRQLAL LHFRNTITLS VKLEDALARS HACVPPAIAQ MLLVLQGVHE SRGVTEYLR LETLIQKVVV PYLGTYGLYS NEGPCTHSCI LEKRFLRRSR SGDILAKNPV VRSKSYNTPL LNPVAEHEAE GTASGGTSIR RHSVSEMTSC PEPQGFVDTP DQGPGTFRS SPSPHSGPCP SRLYPPAHSP EQGPDHGSPP TSSPETLVDQ ILESVDSDSE GIFIDFGRGS RSSVSDFEAA GGRPSVV
<b>Source</b>	Baculovirus
<b>Target Names</b>	Prr5
<b>Protein Names</b>	Recommended name: Proline-rich protein 5 Alternative name(s): Protein observed with Rictor-1 Short name= Protor-1
<b>Expression Region</b>	1-387
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a protein with a proline-rich domain. This gene is located in a region of chromosome 22 reported to contain a tumor suppressor gene that may be involved in breast and colorectal tumorigenesis. Rare read-through transcripts, containing exons from the ARHGAP8 gene which is located immediately downstream, led to the original description of PRR5 and ARHGAP8 as a single gene. Alternative splicing and the use of alternative promoters results in transcripts encoding different isoforms.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.