



# Recombinant Rat Protein FAM98A (Fam98a)

<b>Product Code</b>	CSB-EP688878RA-B
<b>Abbreviation</b>	Fam98a
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5FWT1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MECDLMETDI LESLEDLGYK GPLLDDGALL QAVSAGAASP EFTKLCAWLV SELRVLCKLE ENVQATNSPS EAEEFQLEVS GLLGEMNCPY PCLTSGDVTK RLLVQKNCLL LLTYLISELE AARMLCVNAP PKKAQEGGGS EVFQELKGIC IALGMSKPPA NITMFQFFSG IEKCLKETLA KVPPNHVGKP LLKKSMPAH WEKIEAINQA IANEYEVRRK LLIKRLDVTV QSGWSDRAK SQTEKLAKVY QPKRSVLSPK GNISIAHLLA ARQDLSKILR TSSGSIREKT ACAINKVLMG RVPDRGGRPN EIEPPPEMP PWQKRQDGPQ QQAGGRGGGR GGYEHSSYGG RGGHEQGGRG GRGGYDHGGR GGGRGNKHQG GWTGGSASG GGYQDGGYRD PGFQPGGYHG GHSSGYQGGG YGGFQTSSYT GSGYQGGGYQ QDNRYQDGGH HGERGSGRGG RGGRGGRGGR GSQGGGWGGR GSQTYHQGGQ FEQHFQHGYY QYSHSGFGQG RHYTS
<b>Source</b>	E.coli
<b>Target Names</b>	Fam98a
<b>Protein Names</b>	Recommended name: Protein FAM98A
<b>Expression Region</b>	1-515
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.