



Recombinant Rat Protein FAM98A (Fam98a)

Product Code	CSB-EP688878RA
Abbreviation	Fam98a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5FWT1
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MECDLMETDI LESLEDLGYK GPLLDDGALL QAVSAGAASP EFTKLCAWLV SELRVLCKLE ENVQATNSPS EAEEFQLEVS GLLGEMNCPY PCLTSGDVTK RLLVQKNCLL LLTYLISELE AARMLCVNAP PKKAQEGGGS EVFQELKGIC IALGMSKPPA NITMFQFFSG IEKCLKETLA KVPPNHVGKP LLKKSMPAH WEKIEAINQA IANEYEVRRK LLIKRLDVTV QSGWSDRAK SQTEKLAKVY QPKRSVLSPK GNISIAHLLA ARQDLKILR TSSGSIREKT ACAINKVLMG RVPDRGGRPN EIEPPPEMP PWQKRQDGPQ QQAGGRGGGR GGYEHSSYGG RGGHEQGGGR GRGGYDHGGR GGGRGNKHQG GWTDDGSASG GGYQDGGYRD PGFQPGGYHG GHSSGYQGGG YGGFQTSSYT GSGYQGGGYQ QDNRYQDGGH HGERGSGRGG RGGRGGRGGR GSQGGGWGGR GSQTYHQGGQ FEQHFQHGGY QYSHSGFGQG RHYTS
Source	E.coli
Target Names	Fam98a
Protein Names	Recommended name: Protein FAM98A
Expression Region	1-515
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.