



Recombinant Staphylococcus epidermidis Lipoyl synthase (lipA)

Product Code	CSB-MP689179SAAA
Abbreviation	lipA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5HQN7
Product Type	Recombinant Protein
Immunogen Species	Staphylococcus epidermidis (strain ATCC 35984 / RP62A)
Purity	>85% (SDS-PAGE)
Sequence	MATRNEEILR KPDWLKIKLN TNDNYTGLKK MMREKNLHTV CEEAKCPNIH ECWGARRTAT FMILGAVCTR ACRFCVKTG LPNELDLNEP ERVAESVELM NLKHVVITAV ARDDLQDQGS NVYAETVRKV RERNPFTTIE ILPSDMGGDY EALETLMASR PDILNHNIEV VRRLLTPRVRA RATYDRTLQF LRRSKELQPD IPTKSSLMVG LGETMEEIYE TMDDL RANDV DILTIGQYLQ PSRKHLKVEK YYTPLEFGKM RKIAMEKGFK HCQAGPLVRS SYHADEQVNE AAKEKQRQGE EQLN
Source	Mammalian cell
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.