



Recombinant Human Leucine-rich repeat-containing protein 73 (LRRC73)

Product Code	CSB-BP689338HU
Abbreviation	LRRC73
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JTD7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLPSSIQISG EPLSGAEVRD ICRGLRDNAV RLLSLRGCR L CDRDFGRICR ALAGATSLAQ LNLNLGVVSS PSRIKQLAEA LRTNRSIQSL FLHGSPLTDA GLALLNPALA LHPALVALDL GDCMLGDEAI NLICGLLPPD GAKSGLKELT LSANPGITPK GWSRLAIAVA HSSQVRVLNL DYNPLGDHVA GMLAVAVASS RTLEVL DLEG TGLTNQSAQT LLDMVENYPT ALRSLVLAEN SISPELQQQI CDLLSEGE EEE EEVAGGAGDT QEWERGREPA AHQRGSSSWM CPSDPSSQMV LMTSGLGDSL LAETEM
Source	Baculovirus
Target Names	LRRC73
Protein Names	Recommended name: Leucine-rich repeat-containing protein 73
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.