



Recombinant Human T-complex-associated testis-expressed protein 1 (TCTE1)

Product Code	CSB-BP689339HU
Abbreviation	TCTE1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JU00
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MQDVTTSAL LDPSHSSVST QDNSSTGGHT SSTSPQLSKP SITPVPKSR NPHPRANIRR MRRIIAEDPE WSLAIVPLLT ELCIQHIIRN FQKNPILKQM LPEHQKQVLN HLSPDLPLAV TANLIDSENY WLRCCMHRWP VCHVAHHGGS WKRMFFERHL ENLLKHFIPG TTDPAVILD LPLCRNYVRR VHVDQFLPPV QLPAQLRPGD QSDSGSEGEM EEPTVDHYQL GDLVAGLSHL EELDLYDVK DCGMNFEWNL FLFTYRDCLS LAAAIKACHT LKIFKLTRSK VDDDKARIII RSLLDHPVLE ELDLSQNLIG DRGARGAAKL LSHSRLRVLN LANNQVRAPG AQSLAHALAH NTNLI SLNLR LNCIEDEGGQ ALAHALQTNK CLTTLHLGGN ELSEPTATLL SQVLAIN TTL TSINLSCNHI GLDGGKQLLE GMSDNKTLLE FDLRLSDVAQ ESEYLIGQAL YANREARQR ALNPSHF MST ITANGPEN SV G
Source	Baculovirus
Target Names	TCTE1
Protein Names	Recommended name: T-complex-associated testis-expressed protein 1 Short name= Tcte-1
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.