



Recombinant Flap endonuclease 1 (FEN1)

Product Code	CSB-EP689812TCK-B
Abbreviation	FEN1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4N3S6
Product Type	Recombinant Protein
Immunogen Species	Theileria parva (East coast fever infection agent)
Purity	≥85% (SDS-PAGE)
Sequence	MGIKGLIPFL SEKVPSSISE LSLESLSGES LAIDASAALY QFTIAIRDSS YFSSLVNSKG ESTSHIYGLM NRCSKFLEYG IKPVFVFDSK PPELTKTKLE KRRQQREEAN ASLKKAISEG DKESVKKLVG RTVKVSKEMN ESAKLLRLM GVPVIEALEE AEAQCAYLVT KNLCRFVASE DTDTLVFGGA FLLRNVASSS SKKILKVDLQ KVLDGLEFNF DQFIDFCILC GCDYCDTLEG VGPKTAYSLV KKYQNLEEIV NFKGGDYDDF KEAKEYFLSP KVNEYDENS V KLGTIDPEGL TEFLVQENNF SKERVEKFIE KLLKFKTKKI QTSLLSFLTA PQHNNKAKSS NKRPREP KAL DCKPNTY GSG KGTVNVKKES STIKKDEIDL TAEDLFCEPS NSDNEEDDRG RVDKNEDLFK KSENETNEIK QNQFKSKSDE KKKQEEHHTI EFNFNHRKVI LIDDDDDDEVV TTANSEKMNC GKCLIIYFLQ KLIK
Source	E.coli
Target Names	FEN1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.