



Recombinant Arabidopsis thaliana Agamous-like MADS-box protein AGL61 (AGL61)

Product Code	CSB-YP689852DOA
Abbreviation	AGL61
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4PSU4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MYVTKYKNKL PTLQVDLIIM PHQTQACIYK QTLSLHQIFQ QRKATKTLHT KKTMMSSKKKE SIGRQKIPMV KIKKESHRQV TFSKRRAGLF KKASELCTLC GAEIGIIVFS PAKKPFSGH PSVESVLDRY VSRNNMSLAQ SQQLQGSPAA SCELNMQQLTH ILSEVEEEKK KGQAMEEMRK ESVRRSMINW WEKPVEEMNM VQLQEMKYAL EELRKT VVTN MASFNEAKDD VFGFLDNKVT VPPYVNMPSG PSNIYNFANG NGCF
Source	Yeast
Target Names	AGL61
Protein Names	Recommended name: Agamous-like MADS-box protein AGL61 Alternative name(s): Protein DIANA
Expression Region	1-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.