



# Recombinant Mouse Lymphoid enhancer-binding factor 1 (Lef1)

<b>Product Code</b>	CSB-BP012856MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P27782
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPQLSGGGGG GDPELCATDE MIPFKDEGDP QKEKIFAEIS HPEEEGDLAD IKSSLVNESE IIPASNGHEV VRQAPSSQEP YHDKAREHPD EGKHPDGGGLY NKGPSYSSYS GYIMPMNMS DPYMSNGSLS PPIPRTSNKV PVVQPSHAVH PLTPLITYSD EHFSPGSHPS HIPSDVNSKQ GMSRHPPAPE IPTFYPLSPG GVGQITPIIG WQGQPVPYIPIT GGFRQPYPSS LSGDTSMSRF SHHMIPGPPG PHTTGIPHPA IVTPQVKQEH PHTDSDLMHV KPQHEQRKEQ EPKRPHIKKP LNAFMLYMK E MRANVVAECT LKESAAINQI LGRRWHALSR EEQAKYYELA RKERQLHMQL YPGWSARDNY GKKKKRKREK LQUESTSGTGP RMTAAYI
<b>Source</b>	Baculovirus
<b>Target Names</b>	Lef1
<b>Protein Names</b>	Recommended name: Lymphoid enhancer-binding factor 1 Short name= LEF-1
<b>Expression Region</b>	1-397
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene encodes a transcription factor belonging to a family of proteins that share homology with the high mobility group protein-1. This protein can bind to a functionally important site in the T-cell receptor-alpha enhancer, thereby conferring maximal enhancer activity. This transcription factor is involved in the Wnt signaling pathway, and it may function in hair cell differentiation and follicle morphogenesis. Mutations in this gene have been found in somatic sebaceous tumors. This gene has also been linked to other cancers, including androgen-independent prostate cancer. Alternative splicing results in multiple transcript variants.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.