



Recombinant *Macaca fascicularis* F-box only protein 25 (FBXO25)

Product Code	CSB-EP689937MOV-B
Abbreviation	FBXO25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4R372
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Purity	>85% (SDS-PAGE)
Sequence	MPFLGQDWRS PGWSWIKTED GWKRCESCSQ KLERENNHCHN ISHSIILNSE DGEIFNNEEH EYASKKRKGD HFRNDTNTQS FYREKWIYVH KESTKERHGY CTLGAEAFNRL DFSSAIQDIR RFNYVVKLLQ LIAKSQTLTSL SGVAQKNYFN ILDKIVQKVL DDHHNPRLIK DLLQDLSSTL CILIRGVGKS VLVGNINIWI CRLETILAWQ QQLQDLQMTK QVNNGLTSLD LPLHMLNNIL YRFSDGWDII TLGQVTPTLY MLSEDRQLWK KLCQYHFAEK QFCRHLILSE KGHIEWKLMY FALQKHYPK EQYGDTLHFC RHCSILFWKD SGHPCTAADP DSCFTPVSPQ HFIDLFKF
Source	<i>E. coli</i>
Target Names	FBXO25
Protein Names	Recommended name: F-box only protein 25
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.