



Recombinant Human T-cell leukemia virus 3 Protein Tax-3 (tax)

Product Code	CSB-MP690009HAAJ
Abbreviation	tax
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4U0X7
Product Type	Recombinant Protein
Immunogen Species	Human T-cell leukemia virus 3 (strain Pyl43) (HTLV-3)
Purity	>85% (SDS-PAGE)
Sequence	MAHFPGFGQS LLYGYPVYVF GDCVQADWCP ISGGLCSARL HRHALLATCP EHQITWDPID GRVVSSALQY LIPRLPSFPT QRTRTLKVL TPPTTATTPK VPPSFFHAVK KHTPFRNCL ELTLGEQLPA MSFPDPGLRP QNVYTIWGCS VVCLYLYQLS PPMTWPLIPH VIFCHPEQLG AFLTRVPTKR LEELLYKIFL STGAIILPE NCFPTTLFQP TRAPAIQAPW HTGLLPCQKE IVTPGLIWTF TDGSPMISGP CPKEGQPSLV VQSSTFIFQQ FQTKASHPAF LLSHKLIQYS SFHSLHLLFE EYSTVPFSL FNEKGANVSD DEPRGGPQPP TGGQIAESSV
Source	Mammalian cell
Target Names	tax
Protein Names	Recommended name: Protein Tax-3 Alternative name(s): Trans-activating transcriptional regulatory protein of HTLV-3
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.