



Recombinant Danio rerio Poly (U)-specific endoribonuclease-B (endoub)

Product Code	CSB-EP690355DIL
Abbreviation	endoub
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q503V9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MTEGDRELSA LIQELWDNDT NRLRPGIDYR ISLQ GKAGDL AGVPDDSDGA GMPLFSFVDE NIFKKETFLA FISLLDNYES DTGEPEIVTP EEEMENHRFL DSVIKTPTMK IAHKYLVEKR LSPEDTTGFK QQLYKIWFEL YARKGSSKPD SSGFEHV FVG ETRGGHTVIG FHNWIQLYLQ EKLGHIDYKG YSVNANSPQP DENKHM LALQ FSWKNGIKPK GSIFIGVSPE FEFSLYTLCF LMSPNERVKV SFNLYDVEIV CHHYNQKHIG TYPVPVKYL NV
Source	E.coli
Target Names	endoub
Protein Names	Recommended name: Poly(U)-specific endoribonuclease-B EC= 3.1.-.- Alternative name(s): Protein endoU-B Uridylate-specific endoribonuclease-B
Expression Region	1-282
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.