



Recombinant Protein *klaB* (*klaB*)

Product Code	CSB-MP690419ENL
Abbreviation	<i>klaB</i>
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q52328
Product Type	Recombinant Protein
Immunogen Species	<i>Escherichia coli</i>
Purity	>85% (SDS-PAGE)
Sequence	MNALKTTHDA KAPIVAFDMT PATLRELGLQ ESDVPEVHAV AQRIEVGSPQ TVAEFGRDVA EHTSRYADSL LDQVRNSDL D EAGEKLTQVV AKARSLNVGP LSDNRSRLPL IGPLIDFRV RSTGFMARFD TTREQIEHLV SEVQTTQQGI AQRNASLDEM FAAVREEHRL LGVHIAAGKV RLAELREQAE GLRGNVGNPD GRVQELADLD AMVANLDKRI GDLIALQHSA MQSLPTIRMI QANNQMLV DK FHTIREITVP AWKRQFMLAL SLNEQKNAVE LATAIDDTN DLMKRNAALL HRTSVETAKE NQRLVIDVDT LKQVQTTLIK TVEDVIRIQQ EGVQKRKDAE KQIAAMRGDL QAKLTRQPVR ELAQQESV
Source	Mammalian cell
Target Names	<i>klaB</i>
Protein Names	Recommended name: Protein <i>klaB</i> Alternative name(s): Protein <i>telA</i>
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.