



Recombinant *Oryza sativa* subsp. japonica Myb-related protein MYBAS1 (MYBAS1)

Product Code	CSB-YP690483OFG
Abbreviation	MYBAS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53NK6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MVTVREEMRK GPWTEQEDLQ LVCTVRLFGD RRWDFVAKVS GLNRTGKSCR LRWVNYLHPG LKHGRMSPKE EHLIELHAR WGNRWSRIAR RLPGRTDNEI KNYWRTHMRK KAQERRGDMS PSSSSSSLVY QSCLLDTVPI ISMDGGDIHD DRSCMARVLK STQSVMDGYT MDQIWKEIEA PGAPSLLGID EGKDKACSNL PCPLLTSTMS DYSCPEVFWK IDNEETRMLA TQSGYGK
Source	Yeast
Target Names	MYBAS1
Protein Names	Recommended name: Myb-related protein MYBAS1
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.