



Recombinant Dictyostelium discoideum Ubiquitin-like-conjugating enzyme ATG10 (atg10)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-MP690854DKK |
| Abbreviation | atg10 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q55EL2 |
| Product Type | Recombinant Protein |
| Immunogen Species | Dictyostelium discoideum (Slime mold) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MLTSKDFRDQ AINLIKKWNN IIDEIPWQWN QINELNNE SK GYFTTKRYHK INNNNNNNNN NNIENKNNNN IENFEEIKET IDDSSTTIK SNNNNNENNI IIFQFDIYS KSYQVPVLYL NGFSSFSSP LSWNEIWN NL PLSNLDKNQQ STIPYITQVE HPILGNPCYQ LHPCETDNLM KLILLKEKDY NDNNDKKEYF KDYLLSWLS IIGPMVNIKI PFDLLKNNNI |
| Source | Mammalian cell |
| Target Names | atg10 |
| Protein Names | Recommended name: Ubiquitin-like-conjugating enzyme ATG10 EC= 6.3.2.- Alternative name(s): Autophagy-related protein 10 |
| Expression Region | 1-230 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |