



Recombinant Rat General transcription factor IIH subunit 3 (Gtf2h3)

Product Code	CSB-MP690889RA
Abbreviation	Gtf2h3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q561R7
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MAADEDELNL LVIIVDTNPI WWGKQALKES QFTLSKCMDA VMVLANAHLF MNRSNQLAVI ASHIQESRFL YPGKNGRLGD FFGDPGNALP DCNPSGSKDG KYELLTAANE VIAEEIKDLM TKSDIKGQHT ETLLAGSLAK ALCYIHRASK AVKDNQEMKS RILVIKAAED SALQYMNFMN VIFAAQKQNI LIDACVLDS SGLLQQACDI TGGLYLKVPQ MPSLLQYLLW VFLPDQDQRS QLILPPPIHV DYRAACFCHR SLIEIGYVCS VCLSIFCNFS PICTTCETAF KISLPPVLKA KKKKQKVSL
Source	Mammalian cell
Target Names	Gtf2h3
Protein Names	Recommended name: General transcription factor IIH subunit 3 Alternative name(s): General transcription factor IIH polypeptide 3
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.