



Recombinant Human KHDC3-like protein (KHDC3L)

Product Code	CSB-MP691268HU
Abbreviation	ECAT1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q587J8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDAPRRFPTL VQLMQPKAMP VEVLGHLPKR FSWFHSEFLK NPKVVRLEVV LVEKIFGRGG ERIPHVQGMS QILIHVNRLD PNGEAEILVF GRPSYQEDTI KMIMNLADYH RQLQAKGSGK ALAQDVATQK AETQRSSIEV REAGTQRSVE VREAGTQRSV EVQEVGTQGS PVEVQEAGTQ QSLQAANKSG TQRSPEAASK AVTQRFREDA RDPVTRL
Source	Mammalian cell
Target Names	KHDC3L
Protein Names	Recommended name: ES cell-associated transcript 1 protein
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.