



Recombinant Human RILP-like protein 1 (RILPL1)

Product Code	CSB-MP691645HU
Abbreviation	RILPL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5EBL4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEEERGSALA AESALEKNVA ELTVMDVYDI ASLVGHEFER VIDQHGCEAI ARLMPKVVRV LEILEVLVSR HHVAPELDEL RLELDRLRLE RMDRIEKERK HQKELELVED VWRGEAQDLL SQIAQLQEEN KQLMTNLSHK DVNFSEEEFQ KHEGMSERER QVMKKLKEVV DKQRDEIRAK DRELGLKNED VEALQQQQTR LMKINHDLRH RVTVVEAQGK ALIEQKVELE ADLQTKEQEM GSLRAELGKL RERLQGEHSQ NGEEEPETEP VGEESISDAE KVAMDLKDPN RPRFTLQELR DVLHERNELK SKVFLLEEL AYYKSEEMEE ENRIPQPPPI AHPRTSPQPE SGIKRLFSFF SRDKKRLANT QRNVHIQESF GQWANTHRDD GYTEQGQEAL QHL
Source	Mammalian cell
Target Names	RILPL1
Protein Names	Recommended name: RILP-like protein 1 Alternative name(s): Rab-interacting lysosomal-like protein 1
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.