



Recombinant Human Tumor necrosis factor alpha-induced protein 8-like protein 3 (TNFAIP8L3)

Product Code	CSB-BP691848HU
Abbreviation	TNFAIP8L3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5GJ75
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGKPRQNPST LVSTLCEAEP KGKLWVNGYA GTQGTRDATL QTRLIPLSFH LQRGKGLAAP LSALSAPRLP ERPADGRVAV DAQPAARSMD SDSGEQSEGE PVTAAGPDVF SSKSLALQAQ KKILSKIASK TVANMLIDDT SSEIFDELYK VTKEHTHNKK EAHKIMKDLI KVAIKIGILY RNNQFSQEEL VIVEKFRKKL NQTAMTIVSF YEVEYTFDRN VLSNLLHECK DLVHELVRQH LTPRTHGRIN HVFNFHADVE FLSTLYSLDG DCRPNLKRIC EGINKLLDEK VL
Source	Baculovirus
Target Names	TNFAIP8L3
Protein Names	Recommended name: Tumor necrosis factor alpha-induced protein 8-like protein 3 Short name= TNF alpha-induced protein 8-like protein 3 Short name= TNFAIP8-like protein 3
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.