



Recombinant Pan troglodytes Fibroblast growth factor 2 (FGF2)

Product Code	CSB-EP692217EQV
Abbreviation	FGF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5IS69
Product Type	Recombinant Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Purity	≥85% (SDS-PAGE)
Sequence	PALPEDGG SGAFPPGHFK DPKRLYCKNG GFFLRIHPDG RVDGVREKSD PHIKLQLQAE ERGVVSIKGV CANRYLAMKE DGRLLASKCV TDECFFFERL ESNNYNTYRS RKYTSWYVAL KRTGQYKLGS KTGPGQKAIL FLPMSAKS
Source	E.coli
Target Names	FGF2
Protein Names	Recommended name: Fibroblast growth factor 2 Short name= FGF-2 Alternative name(s): Basic fibroblast growth factor Short name= bFGF Heparin-binding growth factor 2 Short name= HBGF-2
Expression Region	143-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.