



Recombinant *Oryza sativa* subsp. japonica Zinc finger A20 and AN1 domain-containing stress-associated protein 3 (SAP3)

Product Code	CSB-EP692280OFG
Abbreviation	SAP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JN07
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGQQVQHESR INVGEATHVS KAEMGANTMF AT SRLNSNNK VGPELAYLSG VASSASDSST AAPSPCYLCH KPAALHVFG LAGRYVFGSVK REAYLSQEGP RSGRTPNRIA ESLPVRVVND FGLRLRVVTN QGPIKPRPPR PIDAIVFASI ETRNRLRGFD RSFCCSAPPE TYVFLPRARE TIVLRANI IK MSSEQQASAG QPVLCASGCG FYGNPATLDM CSVCYRQHCL LNGATMATGP SSSVAAASAA TVATGAVTSD SCSVPSAEVN GAAFSSKNNP EPATVVEKKA PANRCASCKK KVGLLGFACR CGATYCGTHR YPEKHACGFD FKGASRDAIA RANPLIKGEK LTNKI
Source	E.coli
Target Names	SAP3
Protein Names	Recommended name: Zinc finger A20 and AN1 domain-containing stress-associated protein 3 Short name= OsSAP3
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.