



# Recombinant Human Long-chain specific acyl-CoA dehydrogenase, mitochondrial (ACADL)

<b>Product Code</b>	CSB-EP001125HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P28330
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GGEERLETPS AKKLTDIGIR RIFSPEHDIF RKSVRKFFQE EVIPHHSEWE KAGEVSRE VW EKAGKQGLLG VNIAEHLGGI GGDLYSAAIV WEEQAYSNC S GPGFSIHSGI VMSYITNHGS EEQIKHFIPQ MTAGKCIGAI AMTEPGAGSD LQGIKTNAKK DGSDWILNGS KVFISNGSLS DVVIVVAVTN HEAPSPA HGI SLFLVENG MK GFIKGRKLHK MGLKAQDTAE LFFEDIRLPA SALLGEENKG FYYIMKELPQ ERLLIADVAI SASEFMFEET RNYVKQRKAF GKTVAHLQTV QHKLAELKTH ICVTRAFVDN CLQLHEAKRL DSATACMAKY WASELQNSVA YDCVQLHGGW GYMWEYPIAK AYVDARVQPI YGGTNEIMKE LIAREIVFDK
<b>Source</b>	E.coli
<b>Target Names</b>	ACADL
<b>Protein Names</b>	Recommended name: Long-chain specific acyl-CoA dehydrogenase, mitochondrial Short name= LCAD EC= 1.3.99.13
<b>Expression Region</b>	31-430
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein belongs to the acyl-CoA dehydrogenase family, which is a family of mitochondrial flavoenzymes involved in fatty acid and branched chain amino-acid metabolism. This protein is one of the four enzymes that catalyze the initial step of mitochondrial beta-oxidation of straight-chain fatty acid. Defects in this gene are the cause of long-chain acyl-CoA dehydrogenase (LCAD) deficiency, leading to nonketotic hypoglycemia.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.