



Recombinant Danio rerio UPF0762 protein C6orf58 homolog (zgc:77778)

Product Code	CSB-BP700122DIL
Abbreviation	zgc:77778
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4QRF7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	A VVTENGLPIQ WRKAPSDL SH LPISNTGVQV NPWVYSQRMT MLKMOVINATN AYMSSMGPG E QENPLWSLPL QLGWKLKSGR LADPTPDSSS TCGSEASDPV CISPLSFFGC VNYLSVLPF LAAVETGVVS SGGHQMLIQV PAEVAQDYCS SYSDCSTKHP NTMAKWHLFF QSLRQVSQSK ESDFNKKDSI LGLMWAAEEE SIQTASGACT ERQKLYSSPE VSFQQSSVNL GAFISAAHYH ASIERSGKFM SHLPSRVLQE ADSPPIADL STEENHALYM LSWMNSINQL LGGSLVDLWR KAMCSAQARE KGQAFLDDLI LDPKFP GSSL WSIICVMSTS C
Source	Baculovirus
Target Names	leg1b
Protein Names	Recommended name: UPF0762 protein C6orf58 homolog
Expression Region	20-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.